

Flexibles-Zystoskop

Flexible-Cystoscope



21.000 PIXEL

EXCLUSIVE



Flexibles-Zystoskop Flexible-Cystoscope



Technische Daten Technical Data

Artikel Nr.	/ Item no.	FCS5523.01	FCS5523.11
Einführschlauch	/ Insertion tube	Ø 5.5 mm	Ø 5.5 mm
Arbeitskanal	/ Working channel	Ø 2.3 mm	Ø 2.3 mm
Nutzlänge	/ Working length	380 mm	380 mm
Blickwinkel	/ Field of view	85°	85°
Abwinklung auf/ab	/ Angulation up/down	210°/120°	210°/120°
Bildpunkte	/ Pixel	21.000	21.000
Tiefenschärfe	/ Field of depth	3 - 50 mm	3 - 50 mm
Absaugventil	/ Suction valve	✓	-

Zubehör / accessories	In der Lieferung enthalten / including delivery	Artikel Nr. / Item no.
Dichtigkeitstester / Leak tester	✓	200.00040.02
ETO-Kappe / ETO-Cap	✓	200.02421.00
Absaugventil / Suction valve	✓	370.00249.00
Reinigungsbürste / Cleaning brush	-	VTL/EDB-120-03
Spüladapter / Flushing adapter (FCS5523.01)	✓	200.04203.00
Biopsieventil / Biopsy valve (FCS5523.01)	✓	VTL/BP02S
Hahnbrücke mit Luer-Anschluß und zwei Hähnen / Bridge with Luer-Lock, and two stop cocks (FCS5523.11)	-	220.00003.00
Justierbarer Einweg Y-Spül, Instrumenten-Ansatz / Adjustable Disposable Y-Flush, Instrumentation Valve (FCS5523.11)	✓	200.04105.ce (Einweg / Single use)

Bild / Picture	Artikelnummer / Item No.	Beschreibung / Description
	AFF204-216-060	Flexible Fass-/Fremdkörperzange, Alligatormaul, 600 mm, 5 Charr. Flexible grasping/foreign body forceps, serrated, 600 mm, 5 Fr.
	AFF203-216-060	Flexible Probeexcisions- und Fasszange, gezahnt, 600 mm, 5 Charr. Flexible biopsy and grasping forceps, serrated, 600 mm, 5 Fr.