

## Zystoskope Cystoscopes



# Zystoskope

## Cystoscopes



## Technische Daten

### Technical Data

**Zystoskop** Ø 2.8 mm / **Cystoscope** Ø 2.8 mm

Artikel Nr. / Item no.	US28A.01	US28F.01	US28B.01
Durchmesser / Diameter	Ø 2.8 mm	Ø 2.8 mm	Ø 2.8 mm
Kompatibel / Compatible	Storz®	Storz®	Storz®
Nutzlänge / Working length	302,5 mm	302,5 mm	302,5 mm
Blickrichtung / Direction of view	0°	12°	30°

**Zystoskop** Ø 4.0 mm / **Cystoscope** Ø 4.0 mm

Artikel Nr. / Item no.	US40A.01	US40B.01
Durchmesser / Diameter	Ø 4.0 mm	Ø 4.0 mm
Kompatibel / Compatible	Storz®	Storz®
Nutzlänge / Working length	302,5 mm	302,5 mm
Blickrichtung / Direction of view	0°	30°

Artikel Nr. / Item no.	US40C.01	US40F.01
Durchmesser / Diameter	Ø 4.0 mm	Ø 4.0 mm
Kompatibel / Compatible	Storz®	Storz®
Nutzlänge / Working length	302,5 mm	302,5 mm
Blickrichtung / Direction of view	70°	12°