





Operating Surgical
Shadowless Lamp

TOOTOO MEDITECH CO., LTD.



D5/D5

- ★ New petal style three-dimensional design
- ★ Imported German OSRAM lamp beads
- ★ Imported French lens



Optional with original German imported balance arm

★True cold light source

Adopting a new type of LED cold light source, energy-saving, environmentally friendly, and ultra longThe service life can reach over 80000 hours.

★Excellent deep lighting

Provide adjustable color temperatures of 3700K-5000K, close to natural sunlight, truly reflecting the color of human tissues, fully meeting various surgical lighting needs.

D3/D5



D5 Ceiling mounted



D3 Ceiling mounted



TTOL-D5
Mobile stand



D3 Wall mounted



D5 Wall mounted
Petal single lamp wall mounted



TTOL-D3
Mobile stand

★Advanced Control System

Adopting LCD display button control, it can adjust power switch, illumination, color temperature, etc. to meet the requirements of medical staff for different patient surgical brightness.

★Constant color temperature

LED color temperature is constant, color temperature does not decay, soft and not dazzling, very close to natural sunlight..



LED700/500

- ★ Imported from Germany with original square balance arm
- ★ Imported German OSRAM lamp beads
- ★ Imported French lens




Direct plug-in - suitable for rooms below 2.6 meters

★Excellent shadowless effect

The lamp cap adopts the most scientific radian, The multi-point light source design makes the spot illumination more abundant, and when the surgical lamp is partially obstructed, it can also achieve a perfect shadowless effect.

The entire series of shadowless lights can be equipped with a high-definition camera teaching system, which can be installed with internal, external, and panoramic cameras according to customer requirements.



LED700 Ceiling mounted

★Imported German OSRAM lamp beads

★Imported French lens

★Intelligent operation (with memory function)



LED700Wall mounted



LED700Mobile

LED500Ceiling mounted



LED500Wall mounted



LED500Mobile



Operating Lamp

D480/380European style



- ★ Imported from Germany with original square balance arm
- ★ Imported German OSRAM lamp beads
- ★ Imported French lens

D480European style Ceiling mounted



D380European style Ceiling mounted



P/05 D480A European style



D380A European style



Operating Lamp

D480/380American style



- ★ Imported from Germany with original square balance arm
- ★ Imported German OSRAM lamp beads
- ★ Imported French lens

D480American style Ceiling mounted



D380American style Ceiling mounted



D480A American style



D380A American style



D480/D480 American style



D480/D480 European style



D480 American style Wall mounted



D480 European style Wall mounted



D380 American style Wall mounted



D380 European style Wall mounted





D300 Ceiling mounted



D300 Surgical Lamp

- 1.Illuminance(Lux)10,000-80,000Lux
- 2.Color Temp (K):3700k-5000K
- 3.Color Temp(K):90
- 4.Depth of Focus: $\geq 1200\text{mm}$
- 5.Spot Diameter:200mm
- 6.Brightness Adjustment:1%-100%
- 7.LED ServiceLife: $\geq 50,000\text{h}$
- 8.Supply Voltage:AC100-240V 50/60HZ
- 9.Doctor's Head Temperature $< 1\%$



D300Wall mounted

D300A



D200 Surgical Lamp

- 1.Illuminance(Lux) 30000
- 2.Color Temp (K) 3500-5000
- 3.Spot Diameter $\phi 50\text{mm}$
- 4.Lighting depth $\geq 350\text{mm}$
- 5.Safety type Class B
- 6.Bulb Power 10W



D200 Ceiling mounted



D200Wall mounted

D200A (deep irradiation)



D200A (Fixed type)



D200B (Fixed type)





ZF700/500



ZF500Wall Mounted



ZF500Wall mounted direct insertion model



ZF700 Ceiling mounted



ZF500 Ceiling mounted



(1)Integral reflection shadowless lamp: With high luminous efficiency, it can achieve ideal brightness, color temperature, and luminous efficiency.

(2)There are two halogen lamps in the lamp holder, one of which is a backup bulb, which can automatically activate immediately when the main bulb is damaged.



ZF500 Mobile plug-in type



ZF500Mobile